


**NEW! Disinfecting  
Deodorizing Wipes**



**Virucidal\* / Bactericidal  
Tuberculocidal**

**Pseudomonicidal  
Fungicidal**

A nonwoven disposable cloth containing a stable, low pH formulated DISINFECTANT and DEODORANT for use on hard, non-porous surfaces, in...

*Emergency Medical Settings  
Medical Offices  
Dental Offices  
School Buses  
Recovery Areas  
Laboratories  
Commercial Facilities*

*Hospitals  
Day Care Facilities  
Orthopedics  
Veterinary Facilities  
Newborn Nursery  
Anesthesia / X-ray  
Industrial Facilities*

*Nursing Facilities  
Intensive Care Units  
Spas & Health Clubs  
Janitorial Facilities  
Surgery Units  
Schools  
Clinics*

where control of cross-contamination is required

*For use on hard non-porous surfaces only, such as-  
STAINLESS STEEL, FORMICA, GLASS TABLES,  
CARTS, BASKETS, COUNTERS, CABINETS,  
and TELEPHONES.*

**Total! Timely! and Terrific!**



**KILLS:**  
**HERPES SIMPLEX  
VIRUS TYPE 1, TYPE 2,  
HIV-1 (AIDS VIRUS)\*  
INFLUENZA  
MRSA**

**50**  
**Pre-saturated  
Wipes  
(7 x 10 inches)**

# aero CleanSwipe Disinfecting/ Deodorizing/Cleaning Wipes

Virucidal\* / Bactericidal  
Tuberculocidal      Pseudomonicidal  
Fungicidal

THE ACTIVITY OF CLEANSWIPE  
DISINFECTING / DEODORIZING /  
CLEANING WIPES AGAINST THE HIV-1 (AIDS VIRUS):

## CLEANSWIPE

**DISINFECTING/DEODORIZING/CLEANING WIPES** completely inactivate the HIV-1 (AIDS virus) on hard, non-porous surfaces when used at a contact time of 1 minute at room temperature (20-25°C).

## CLEANSWIPE

**DISINFECTING/DEODORIZING/CLEANING WIPES** kill HIV-1 on pre-cleaned environmental surfaces/objects soiled with blood/body fluids, in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immune deficiency virus Type-1 (HIV-1) (associated with AIDS).

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:**

**Personal Protection:** Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks or eye coverings.

**Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of **CLEANSWIPE** DISINFECTING/DEODORIZING/CLEANING WIPES.

**Contact Time:** HIV-1 VIRUS IS INACTIVATED IN 1 MINUTE.

**Disposal of Infectious Materials:** Materials that may contain blood or other body fluids should be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal.

**PREPARATION OF DISPENSER:** Remove cover from canister and discard covering seal. From center of wipes roll, pull up a wipe corner, twist into a point and thread through the X hole in the underside of the lid cover. Pull wipe through about one (1) inch. Replace the lid cover on container, assuring a secure fit. Pull out first wipe and snap at a 90° angle. Remaining wipes feed automatically for next use. Always reseal hinged cap after each use to insure moisture retention.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

**DIRECTIONS FOR USE:** It is a violation of U.S. Federal Law to use this Product in a manner inconsistent with its labeling.

**DISINFECTION:** Thoroughly wet pre-cleaned, hard, non-porous surface with a wipe, keep wet for 5 minutes, and allow to air dry. For use as a general disinfectant wipe on hard, nonporous, inanimate surfaces such as chairs, baskets, counter tops, exam tables, cages, toilets, basins, tubs, bed springs, frames, bedside and office furniture, chart racks and equipment, cabinets, door knobs, office machinery, sinks, mirrors, light fixtures and switches, supply, dressing, and trash carts, telephones, waste cans and laundry hampers, wheel chairs, and walkers, urinals, plastic and other similar hard surfaces. You may use a wipe to pre-clean the surface to be disinfected, but thorough re-wiping with a fresh wipe is required for disinfection. Not to be used on food contact surfaces.

**CLEANSWIPE** Disinfecting/Deodorizing/Cleaning Wipes are effective against the following pathogens in 5 minutes at 20°C/68°F, on hard, non-porous surfaces when used as directed: *Staphylococcus aureus*, *Trichophyton mentagrophytes* (Athlete's foot fungus), *Pseudomonas aeruginosa*, *Mycobacterium tuberculosis*, *Salmonella choleraesuis*, *Methicillin Resistant Staphylococcus aureus* (MRSA), *Escherichia coli* 0157:H 7, *Vancomycin Resistant Enterococcus faecalis* (VRE), *\*Herpes simplex virus 2* (genital Herpes virus), *\*HIV-1 (AIDS virus)\*\** (At an exposure time of 1 min.), *\*Influenza A (A<sub>2</sub> Japan)*, *\*\**See panel for complete HIV-1 (AIDS virus) information.

**DEODORIZING:** Wipe hard, non-porous surface thoroughly to eliminate odors from bathrooms, industrial plant processes, fecal matter, urine, burns, smoke, mildew and mold, decay, stale air, tobacco, sickrooms and musty closets.

**CLEANING:** To clean surfaces, wipe thoroughly with a fresh wipe, which may be used to deodorize at the same time as cleaning. Use in bathrooms, industrial plants, kitchens (non-food use surfaces), sickrooms, musty closets and similar areas. When used to clean a hard non-porous surface prior to disinfection, a second wipe must be used for disinfection.

## ACTIVE INGREDIENTS

N-Alkyl (60% C <sub>14</sub> , 30% C <sub>16</sub> , 5% C <sub>12</sub> , 5% C <sub>18</sub> )	
Dimethyl Benzyl Ammonium Chloride.....	0.12%
N-Alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> )	
Dimethyl Ethyl Ammonium Chloride.....	0.12%
Isopropyl Alcohol.....	41.58%
Inert Ingredients†.....	58.18%
TOTAL.....	100.00%

†Does not include weight of dry wipe.

WA40/MAR/03



ATLANTA 770.968.9ABC / 1.800.795.9ABC / fax 770.968.7281  
DALLAS (FT. Worth) 972.988.9200 / fax 972.641.8170